

## Aspire Oxygen Bottle Holder



- 140 mm diameter stainless steel holder suits sizes Oxygen Bottle sizes C to D
- Adaptable clamp fits most shower commode brands

### Warranty Information

- 1 Year



✓ Contracted DVA product

✓ Contract approved

### Support @ Home

**Function:** Washing, bathing and showering

**Category:** Prescribed

Click [here](#) for more information

### Compliance

AS/NZS 3973:2009

| Code      | Display Name          | Material        |
|-----------|-----------------------|-----------------|
| BTC067200 | Oxygen Bottle Carrier | Stainless Steel |