


## Liko Sabina II EE Lifter



- Conventional patient lifting to a limited extent, lifting in the sitting position
- Powerful lifting motor ensures the smooth raising and lowering of patients
- Base-width adjustment offers greater flexibility of use when lifting from a bed or chair, enhancing patient mobility
- The lower leg support and the foot rest can easily be removed for cleaning
- Adjustable height settings accommodate a variety of different patient sizes
- Sling bars are sold separately to give customer the choice of best one for their needs. We recommend either the Sabina II 350 SlingBar (LSS393270) or Sabina II Comfort Sling Bar (LSS393280)

### Warranty Information

- Frames: 10 Years
- Electrical Parts: 2 Years

 Contract approved

**ndis**

**Item Code:**

05\_123604167\_0105\_1\_2

**Hoist –**

**Composite:**05\_420\_0105\_1\_2

Click [here](#) for more information

**Compliance**

AS/NZS ISO 10535:2011

| Code      | Display Name                       | Total Length | Total Width   | Total Height  | Safe Working Load | Material |
|-----------|------------------------------------|--------------|---------------|---------------|-------------------|----------|
| LSS393250 | Liko Sabina II EE Lifter           | 1060 mm      | 695 - 1115 mm | 900 - 1000 mm | 200 kg            | Steel    |
| LSS393320 | Liko Sabina Seat Strap SlingBar    | -            | -             | -             | -                 | -        |
| LSS393330 | Liko Sabina II Seat Strap SlingBar | -            | -             | -             | -                 | -        |
| LSS393340 | Liko Sabina Seat Strap             | -            | -             | -             | -                 | -        |
| LSS393350 | Liko Sabina Calf Strap             | -            | -             | -             | -                 | -        |
| LSS393270 | Liko Sabina II 350 SlingBar        | -            | -             | -             | -                 | -        |

| Code      | Display Name                    | Total Length | Total Width | Total Height | Safe Working Load | Material |
|-----------|---------------------------------|--------------|-------------|--------------|-------------------|----------|
| LSS393280 | Liko Sabina II Comfort SlingBar | -            | -           | -            | -                 | -        |
| LSS393290 | Liko Sabina Side Support        | -            | -           | -            | -                 | -        |
| LSS393300 | Liko Sabina II Side Support     | -            | -           | -            | -                 | -        |
| LSS393310 | Liko Sabina II Heel Support     | -            | -           | -            | -                 | -        |