

Liko Original Highback Sling




- Can be used to lift most patients
- Good choice for lifting from a horizontal position
- Provides support for the entire back and head
- Enables a semi-reclined sitting posture

Warranty Information

- 1 Year



 Contract approved

ndis

Item Code:

05_123603164_0105_1_2

Hoist –

Composite:05_420_0105_1_2

Click [here](#) for more information

Compliance

AS/NZS ISO 10535:2011

| Code | Display Name | Safe Working Load | Material |
|-----------|----------------------------------------------------|-------------------|---------------|
| LSS367854 | Original High Back MOD 200/210 Net Poly - MEDIUM | 200 kg | Net Polyester |
| LSS367856 | Original High Back MOD 200/210 - Net Poly - XLARGE | 500 kg | Net Polyester |