

Wheelchair Transport Summary



The table below summarises the transport testing and configuration for Aidacare wheelchairs.
 If the chair you have (or are considering purchasing) is not on this list then it is not crash-tested which means it is not to be used for transport with user in it.
 For complex rehab chairs, please consult your equipment specialist.

Code	Item Description & Image	Image	Crash Tested	Transport	Headrest Code
ASPIRE EVOKE2 - ADD IMAGE					
MWS449560	WChair Folding - Lightweight - Aspire Evoke2 JUNIOR - 300mm - 150kg - White		Yes	Yes - With headrest	WCA848479
MWS449570	WChair Folding - Lightweight - Aspire Evoke2 Youth 350mm - 150kg - Silver		Yes	Yes - With headrest	WCA848479
MWS449500	WChair Folding - Lightweight - Aspire Evoke2 - 350mm - 150kg - Silver		Yes	Yes - With headrest	WCA848479
MWS449510	WChair Folding - Lightweight - Aspire Evoke2 - 400mm - 150kg - Silver ^		Yes	Yes - With headrest	WCA848480
MWS449520	WChair Folding - Lightweight - Aspire Evoke2 - 450mm - 150kg - Silver ^		Yes	Yes - With headrest	WCA848480
MWS449530	WChair Folding - Lightweight - Aspire Evoke2 - 500mm - 150kg - Silver ^		Yes	Yes - With headrest	WCA848480
MWS449540	WChair Folding - Lightweight - Aspire Evoke2 HD - 570mm - 190kg - Silver ^		Yes	Yes - With headrest	WCA848481
MWS449550	WChair Folding - Lightweight - Aspire Evoke2 HD - 610mm - 190kg - Silver ^		Yes	Yes - With headrest	WCA848481
BREEZY BASIX & RUBIX - ADD IMAGE					
MWS449130	Wheelchair Breezy Rubix2 HD 170kg - 20in with 24in Solid Whls BLUE		Yes	Yes - With headrest	WCA858019
MWS449135	Wheelchair Breezy Rubix2 HD 170kg - 22in with 24in Solid Whls BLUE		Yes	Yes - With headrest	WCA858019
MWS449140	Wheelchair Breezy Rubix2 HD 170kg - 24in with 24in Solid Whls BLUE		Yes	Yes - With headrest	WCA858019
REHAB RX - ADD IMAGE					
MWS449630	WChair - Tilt in Space - Aspire Rehab RX JUNIOR - 320mm - 20in Wheels - 180kg - White		Yes	Yes - In standard config	N/A
MWS449710	WChair - Tilt in Space - Aspire Rehab RX - 410mm - 24in Wheels - 180kg - Black ^		Yes	Yes - In standard config	N/A
MWS449720	WChair - Tilt in Space - Aspire Rehab RX - 460mm - 24in Wheels - 180kg - Black ^		Yes	Yes - In standard config	N/A
MWS449740	WChair - Tilt in Space - Aspire Rehab RX - 510mm - 24in Wheels - 180kg - Black ^		Yes	Yes - In standard config	N/A
MWS449750	WChair - Tilt in Space - Aspire Rehab RX2 JUNIOR - 320mm - 20in Wheels - 150kg - White		Yes	Yes - In standard config	N/A
MWS449760	WChair - Tilt in Space - Aspire Rehab RX2 - 410mm - 24in Wheels - 150kg - Black		Yes	Yes - In standard config	N/A
MWS449770	WChair - Tilt in Space - Aspire Rehab RX2 - 460mm - 24in Wheels - 150kg - Black		Yes	Yes - In standard config	N/A
MWS449780	WChair - Tilt in Space - Aspire Rehab RX2 - 510mm - 24in Wheels - 150kg - Black		Yes	Yes - In standard config	N/A